Anti-Canine PD-1(INTERVET 4F12) mAb Cat. No. BME100202



## **PRODUCT INFORMATION**

PD1, PDCD1, CD279, SLEB2, hPD-1, hPD-l, hSLE1
ELISA, Flow Cyt
ELISA 1:5000-10000, Flow Cyt 1:100
Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Chimeric Rabbit
IgG
Canine
PD-1
NP_001301026.1
Anti-Canine PD-1(INTERVET 4F12) mAb
In Stock
Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature. Research grade biosimilar. Not for use in
therapeutic or diagnostic procedures for humans or animals.
Research use only
Unconjugated

Email: info@dimabio.com Website: www.dimabio.com







0.2 µg of Canine PD-1, His tagged protein per well




Figure 1. ELISA plate pre-coated by 2 μg/mL (100 μL/well) Canine PD-1 Protein, His Tag (PME-D100005) can bind Anti-Canine PD-1(INTERVET 4F12) mAb (BME100202) in a linear range of 0.64–80 ng/mL.

## Anti-Canine PD-1(INTERVET 4F12) mAb ELISA

0.2 μg of Canine PD-1, hFc tagged protein per well




Figure 2. ELISA plate pre-coated by 2 μg/mL (100 μL/well) Canine PD-1 Protein, hFc Tag (PME-D100008) can bind Anti-Canine PD-1(INTERVET 4F12) mAb (BME100202) in a linear range of 0.64-80 ng/mL.





## Anti-Canine PD-1(INTERVET 4F12) mAb ELISA

 $0.2~\mu g$  of Canine PD-1, mFc tagged protein per well




Figure 3. ELISA plate pre-coated by 2  $\mu$ g/mL (100  $\mu$ L/well) Canine PD-1 Protein, mFc Tag (PME-D100006) can bind Anti-Canine PD-1(INTERVET 4F12) mAb (BME100202) in a linear range of 0.64–80 ng/mL.




Figure 4. Flow cytometry analysis with 1µg/mL Anti-Canine PD-1(INTERVET 4F12) mAb (BME100202) on Expi293 cells transfected with Canine PD-1 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

Email: info@dimabio.com Website: www.dimabio.com

