

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | LVGN-6051;LVGN6051 |
| Synonyms | ILA;4-1BB;CD137;CDw137;IMD109□TNFRSF9 |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | Whole mAb |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | 4-1BB |
| Uniprot ID | Q07011 |
| Description | Anti-4-1BB(exlinkibart biosimilar□IgG1) mAb |
| Delivery | 3-4 weeks |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

