

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | VIS-649;VIS649;VIS 649 |
| Synonyms | APRIL;CD256;TALL2;ZTNF2;TALL-2;TNLG7B;TRDL-1;UNQ383/PRO715 |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | Whole mAb |
| IgG type | IgG2 |
| Reactivity | Human |
| Target | APRIL |
| Uniprot ID | O75888 |
| Description | Anti-APRIL(sibeprenlimab biosimilar) IgG2) mAb |
| Delivery | 3-4 weeks |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | Unconjugated |

DIMABIO CONFIDENTIAL

