

PRODUCT INFORMATION

Common Name	ADCT-601;ADCT601;ADCT 601;BGB 601;BGB-601;BGB601
Synonyms	AXL;ARK;UFO;AXL3;JTK11;Tyro7
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	AXL
Uniprot ID	P30530
Description	Anti-AXL(mipasetamab biosimilar) IgG1) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	Unconjugated

DIMA BIO PHARMACEUTICAL

