

**PRODUCT INFORMATION**

<b>Common Name</b>	MOR 202;MOR202;MOR-202;MOR-03087;MOR 03087;MOR03087;MOR 3087;MOR-3087;MOR3087;TJ-202;TJ202;TJ 202
<b>Synonyms</b>	CD38;ADPRC1;cADPR1;ADPRC 1;T10;cADPr hydrolase 1
<b>Applications</b>	N/A
<b>Recommended Dilutions</b>	N/A
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Host Species</b>	Whole mAb
<b>IgG type</b>	IgG1
<b>Reactivity</b>	Human
<b>Target</b>	CD38
<b>Uniprot ID</b>	P28907
<b>Description</b>	Anti-CD38(felzartamab biosimilar)IgG1) mAb
<b>Delivery</b>	3-4 weeks
<b>Storage &amp; Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated

DIMA BIO  
CONFIDENTIAL