

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | PF-04236921;PF-4236921;PF 04236921;PF 4236921;PF04236921;PF4236921;TOUR-006;TOUR006;TOUR 006 |
| Synonyms | CDF;HGF;HSF;BSF2;IL-6;BSF-2;IFNB2;IFN-beta-2 |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | Whole mAb |
| IgG type | IgG2 |
| Reactivity | Human |
| Target | IL6 |
| Uniprot ID | P05231 |
| Description | Anti-IL6(pacibekitug biosimilar) IgG2) mAb |
| Delivery | 3-4 weeks |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | Unconjugated |

DIMABIO CONFIDENTIAL

