

PRODUCT INFORMATION

Common Name	MK-0482-Hu52B8
Synonyms	ILT3; LIR5; CD85k
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Humanized
IgG type	IgG4
Reactivity	Human
Target	LILRB4
Uniprot ID	Q8NHJ6
Description	Anti-LILRB4(MK-0482 biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMA BIO CONFIDENTIAL






Figure 1. Flow cytometry analysis with 1µg/mL Anti-LILRB4(MK-0482 biosimilar) mAb (BME100237) on Expi293 cells transfected with Human LILRB4 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

DIMABIO CONFIDENTIAL

