

PRODUCT INFORMATION

Common Name	MNRP1685A
Synonyms	Neuropilin-1;CD304
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Homo sapiens
IgG type	IgG1
Reactivity	Human
Target	NRP1
Uniprot ID	O14786
Description	Anti-NRP1(vesencumab biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 1 µg/mL Anti-NRP1 (vesencumab biosimilar) mAb (BME100119) on Expi293 cells transfected with Human NRP1 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

