

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | ANB 011;ANB011;ANB-011;GSK-4057190;GSK4057190;JEMPERLI;TSR 042;TSR042;TSR-042;WBP 285;WBP285;WBP-285 |
| Synonyms | PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1 |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | Whole mAb |
| IgG type | IgG4 |
| Reactivity | Human |
| Target | PD-1 |
| Uniprot ID | Q15116 |
| Description | Anti-PD-1(dostarlimab biosimilar) IgG4) mAb |
| Delivery | 3-4 weeks |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | Unconjugated |

DIMABIO CONFIDENTIAL

