

**PRODUCT INFORMATION**

|   |   |
|---|---|
| <b>Common Name</b>                      | AB 122;AB-122;AB122;GLS 010;GLS010;GLS010;GLS10;GLS-10;GSL 10;GS-0122;GS0122;GS 0122;GS-122;GS122;GS 122;WBP 3055;WBP3055;WBP-3055;YUTUO  |
| <b>Synonyms</b>                         | PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1  |
| <b>Applications</b>                     | N/A   |
| <b>Recommended Dilutions</b>            | N/A   |
| <b>Formulation &amp; Reconstitution</b> | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.  |
| <b>Host Species</b>                     | Whole mAb   |
| <b>IgG type</b>                         | IgG4  |
| <b>Reactivity</b>                       | Human   |
| <b>Target</b>                           | PD-1  |
| <b>Uniprot ID</b>                       | Q15116  |
| <b>Description</b>                      | Anti-PD-1(zimberelimab biosimilar-IgG4) mAb   |
| <b>Delivery</b>                         | 3-4 weeks   |
| <b>Storage &amp; Shipping</b>           | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| <b>Background</b>                       | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.   |
| <b>Usage</b>                            | Research use only   |
| <b>Conjugate</b>                        | Unconjugated  |

