

PRODUCT INFORMATION

Common Name	AB 122;AB-122;AB122;GLS 010;GLS010;GLS010;GLS10;GLS-10;GSL 10;GS-0122;GS0122;GS 0122;GS-122;GS122;GS 122;WBP 3055;WBP3055;WBP-3055;YUTUO
Synonyms	PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	PD-1
Uniprot ID	Q15116
Description	Anti-PD-1(zimberelimab biosimilar-IgG4) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	Unconjugated

