


**PRODUCT INFORMATION**

<b>Common Name</b>	SEG101, SelG1, crizanlizumab-tmca
<b>Synonyms</b>	P-Selectin;CD62P;SELP;GMP-140
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Host Species</b>	Humanized
<b>IgG type</b>	IgG2
<b>Reactivity</b>	Human
<b>Target</b>	SELP
<b>Uniprot ID</b>	P16109
<b>Description</b>	Anti-SELP(crizanlizumab biosimilar) mAb
<b>Delivery</b>	In Stock
<b>Storage &amp; Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated





**Figure 1.** HEK293 cell line transfected with irrelevant protein (red histogram) and human SELP protein (blue histogram) were surface stained with Anti-SELP(crizanlizumab biosimilar) mAb 1 $\mu$ g/ml followed by PE-conjugated Goat anti-Human IgG secondary antibody.

