

PRODUCT INFORMATION

Common Name	AL-001;AL001;AL 001;AL1;AL-1;AL 1;GSK-4527223;GSK4527223;GSK 4527223
Synonyms	NT3;Gp95;NTR3;LDLCQ6
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	SORT1
Uniprot ID	Q99523
Description	Anti-SORT1(latozinemab biosimilar) IgG1 mAb
Delivery	3-4 weeks
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	Unconjugated

DIMA BIO PHARMACEUTICALS CONFIDENTIAL

