

PRODUCT INFORMATION

Common Name	GEN1029;GEN-1029;Hx-DR5-01;Hx-DR5-05;iDD 004;iDD-004;iDD004
Synonyms	DR5;CD262;KILLER;TRICK2;TRICKB;ZTNFR9;TRAILR2;TRICK2A;TRICK2B;TRAIL-R2;KILLER/DR5
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	TNFRSF10B
Uniprot ID	O14763
Description	Anti-TNFRSF10B(benufutamab biosimilar[IgG1]) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

