

PRODUCT INFORMATION

Target CD46, CD73, CD137, ICAM2, EPCAM, FZD10, GPR33, GUCY2C, HER3
Description The unseeded and seeded tumor sections (Apm sections) from Balb/C nu mice implanted with HT55 cells
Host Cells HT55
Applications IHC
Host mice Balb/C nu
Delivery 1 week
Guarantee 30 months from the date of manufacture
Background These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics.
Storage & Shipping Room temperature
Usage Research use only

DIMABIO CONFIDENTIAL

