

PRODUCT INFORMATION

Target	CD46, CD73, CDH17, CEACAM5, EPCAM, FZD10, GPA33, GUCY2C, HER3
Description	The embedded and sliced tumor sections (4µm section) from Balb/C nu mice implanted with HT55 cells
Host Cells	HT55
Applications	IHC
Host mice	Balb/C nu
Delivery	1 week
Guarantee	36 months from the date of manufacture
Background	These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics.
Storage & Shipping	Room temperature
Usage	Research use only

DIMABIO CONFIDENTIAL

