

PRODUCT INFORMATION

Target	AXL, CD43, EPAM, L6AL51
Description	The embedded and sliced tumor sections (4µm sections) from Balb/C nu mice implanted with PC3 cells.
Host Cells	PC3
Applications	IHC
Host mice	Balb/C nu
Delivery	1 week
Guarantee	30 months from the date of manufacture
Background	These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics.
Storage & Shipping	Room temperature
Usage	Research use only

DIMABIO CONFIDENTIAL

