

PRODUCT INFORMATION

Common Name	BFCR-4350-A, BFCR-4350A, RG-6160, RO-7187797
Synonyms	BXMAS1;CD307;CD307e;FCRH5;IRTA2;PRO820
Applications	ELISA; Flow Cyt
Recommended Dilutions	ELISA 1:5000-10000; Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	FCRL5
Uniprot ID	Q96RD9
Description	Biotinylated Anti-FCRL5(cevostamab biosimilar without anti-CD3 portion) mAb
Delivery	2-3 weeks
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	Biotinylated

