

**PRODUCT INFORMATION**

<b>Common Name</b>	S-78989;S78989;S 78989;VPM 087;VPM087;VPM-087;VPM87;VPM-87;XMA-005.2;XMA-0052;XMA52;XOMA-052;XOMA052;XOMA52;XOMA-52
<b>Synonyms</b>	IL-1;IL1F2;IL1beta;IL1-BETA
<b>Applications</b>	N/A
<b>Recommended Dilutions</b>	N/A
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Host Species</b>	Whole mAb
<b>IgG type</b>	IgG2
<b>Reactivity</b>	Human
<b>Target</b>	IL1B
<b>Uniprot ID</b>	P01584
<b>Description</b>	Biotinylated Anti-IL1B(gevokizumab biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage &amp; Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	Biotinylated

DIMABIO CONFIDENTIAL

