

PRODUCT INFORMATION

Common Name	AEVI-002;AEVI002;AEVI 002;AEVI-2;AEVI2;AEVI 2;AVTX-002;AVTX002;AVTX 002;AVTX-2;AVTX2;AVTX 2;CERC 002;CERC-002;CERC002;CERC 2;CERC-2;CERC2;MDGN-002;MDGN002;MDGN 002;MDGN 2;MDGN-2;MDGN2;SAR 252067;SAR252067;SAR-252067
Conjugate	Biotinylated
Synonyms	LTg;CD258;HVEM;LIGHT;TNFSF14
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG4
Reactivity	Human
Target	LIGHT
Uniprot ID	O43557
Description	Biotinylated Anti-LIGHT(quisovalimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

