

**PRODUCT INFORMATION**

<b>Common Name</b>	MK-4830
<b>Synonyms</b>	LIR-2;ILT-4;MIR-10;CD85d
<b>Applications</b>	ELISA; Flow Cyt
<b>Recommended Dilutions</b>	ELISA 1:5000-10000; Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Host Species</b>	Humanized
<b>IgG type</b>	IgG4
<b>Reactivity</b>	Human
<b>Target</b>	LILRB2
<b>Uniprot ID</b>	Q8N423
<b>Description</b>	Anti-LILRB2 (MK-4830) mAb
<b>Delivery</b>	2-3 weeks
<b>Storage &amp; Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

