

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | BAX069, BAX69 |
| Synonyms | GIF;GLIF;MMIF |
| Applications | ELISA |
| Recommended Dilutions | ELISA 1:5000-10000 |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | Homo sapiens |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | MIF |
| Uniprot ID | P14174 |
| Description | Anti-MIF(imalumab biosimilar) mAb |
| Delivery | 2-3 weeks |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

