

PRODUCT INFORMATION

Clone ID	58A12
Target	TweakR
Synonyms	TNFRSF12A;FGF-inducible 14;FN14;TweakR;CD266
Host Species	Rabbit
Description	Biotinylated Anti-TweakR antibody(58A12), IgG1 Chimeric mAb
Delivery	2-3 weeks
Uniprot ID	Q9NP84
IgG type	Rabbit/Human Fc chimeric IgG1
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1/100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 16% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -88°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	TNFRSF12A (TNF Receptor Superfamily Member 12A) is a Protein Coding gene. Diseases associated with TNFRSF12A include Multiple Sclerosis and Stroke;Ischemic. Among its related pathways are ERK Signaling and CDK-mediated phosphorylation and removal of Cdc6
Usage	Research use only
Conjugate	Biotinylated

