

PRODUCT INFORMATION

| | |
|---|--|
| Target | GP171 |
| Synonyms | H963 |
| Description | Human GP171 full length protein-synthetic nanodisc |
| Delivery | 6~8weeks |
| Uniprot ID | O14626 |
| Expression Host | HEK293 |
| Protein Families | GPCR,Transmembrane,Druggable Genome, |
| Protein Pathways | N/A |
| Molecular Weight | The human full length GP171 protein has a MW of 36.8kDa |
| Formulation & Reconstitution | Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Do not use solvents with pH lower than 6.5 in subsequent experiments. Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Storage & Shipping | |
| Background | Orphan receptor.[UniProtKB/Swiss-Prot Function] |
| Usage | Research use only |

