

PRODUCT INFORMATION

Target	RL3R1
Synonyms	GPCR135, RLN3R1, RXFPR3, SALPR
Description	Human RL3R1 full length protein-synthetic nanodisc
Delivery	6~8weeks
Uniprot ID	Q9NSD7
Expression Host	HEK293
Protein Families	GPCR,Transmembrane,Druggable Genome,
Protein Pathways	GPCRDB Other,Metabolic and Obesity,
Molecular Weight	The human full length RL3R1 protein has a MW of 51.1kDa
Formulation & Reconstitution	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Do not use solvents with pH lower than 6.5 in subsequent experiments.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Receptor for RNL3/relaxin-3. Binding of the ligand inhibit cAMP accumulation.[UniProtKB/Swiss-Prot Function]
Usage	Research use only

