

**PRODUCT INFORMATION**

<b>Target</b>	CCR8, CD45, LGALS1
<b>Description</b>	The embedded and sliced tumor sections (4µm section) from M-NSG mice implanted with Hut78 cells
<b>Host Cells</b>	Hut78
<b>Applications</b>	Identification of immunohistochemical function
<b>Host mice</b>	M-NSG
<b>Delivery</b>	One week
<b>Guarantee</b>	36 months from the date of manufacture
<b>Background</b>	These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics.
<b>Storage &amp; Shipping</b>	Room temperature
<b>Usage</b>	Research use only

DIMABIO CONFIDENTIAL

