

PRODUCT INFORMATION

| | |
|-------------------------------|--|
| Target | CD63 |
| Description | The cryopreserved and sliced tumor sections (4µm sections) from M-NSG mice implanted with T47-D cells. |
| Host Cells | T47-D |
| Applications | IHC |
| Host mice | M-NSG |
| Delivery | 1 week |
| Guarantee | 30 months from the date of manufacture |
| Background | These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics. |
| Storage & Shipping | Room temperature |
| Usage | Research use only |

DIMABIO CONFIDENTIAL

