

PRODUCT INFORMATION

Common Name	LVGN-6051;LVGN6051
Synonyms	ILA;4-1BB;CD137;CDw137;IMD109□TNFRSF9
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	4-1BB
Uniprot ID	Q07011
Description	PE-conjugated Anti-4-1BB(exlinkibart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

