

**PRODUCT INFORMATION**

<b>Common Name</b>	H8, mAb5T4
<b>Conjugate</b>	PE-conjugated
<b>Synonyms</b>	M6P1;WAIF1;TPBG
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	1µl/test
<b>Formulation &amp; Reconstitution</b>	Liquid□PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	Humanized
<b>IgG type</b>	IgG1
<b>Reactivity</b>	Human
<b>Target</b>	5T4
<b>Uniprot ID</b>	Q13641
<b>Description</b>	PE-conjugated Anti-5T4 (H8 biosimilar) mAb
<b>Delivery</b>	In Stock
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>DIMA Disclaimer</b>	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.






Figure 1. Flow cytometry analysis with 1  $\mu$ l/test PE-conjugated Anti-5T4 (H8 biosimilar) mAb on Expi293 cells (Blue histogram) or Jurkat cells (Red histogram).

DIMABIO CONFIDENTIAL

