

PRODUCT INFORMATION

Common Name	DS-6016A;DS6016A;DS 6016A
Synonyms	FOP;ALK2;SKR1;TSRI;ACTRI;ACVR1A;ACVRLK2
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	ACVR1
Uniprot ID	Q04771
Description	PE-conjugated Anti-ACVR1(prafnosbart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

