

PRODUCT INFORMATION

Common Name	MAB-A(Immunogen Inc)
Synonyms	CORD9;MCMP;MDC9;Mltng
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	N/A
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	ADAM9
Uniprot ID	Q13443
Description	PE-conjugated Anti-ADAM9 (biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

