

PRODUCT INFORMATION

Common Name	LY 3002813;LY3002813;LY-3002813
Conjugate	PE-conjugated
Synonyms	AAA;AD1;PN2;ABPP;APPI;CVAP;ABETA;PN-II;preA4;CTFgamma;alpha-sAPP
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	N/A
IgG type	IgG1
Reactivity	Human
Target	APP
Uniprot ID	P05067
Description	PE-conjugated Anti-APP(donanemab biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

DIMABIO CONFIDENTIAL

