

PRODUCT INFORMATION

Common Name	VIS-649;VIS649;VIS 649
Synonyms	APRIL;CD256;TALL2;ZTNF2;TALL-2;TNLG7B;TRDL-1;UNQ383/PRO715
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG2
Reactivity	Human
Target	APRIL
Uniprot ID	O75888
Description	PE-conjugated Anti-APRIL(sibeprenlimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

