

PRODUCT INFORMATION

Common Name	DS-7300;DS7300;DS 7300
Synonyms	B7H3;B7-H3;B7RP-2;4Ig-B7-H3;CD276
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	B7-H3
Uniprot ID	Q5ZPR3
Description	PE-conjugated Anti-B7-H3(ifinatamab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

