

PRODUCT INFORMATION

Common Name	AZD-8205; AZD 8205; AZD8205
Synonyms	B7X;B7H4;B7S1;B7-H4;B7h.5;VCTN1;PRO1291
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	B7-H4
Uniprot ID	Q7Z7D3
Description	PE-conjugated Anti-B7-H4(puxitatug biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

