

PRODUCT INFORMATION

Common Name	SGNB7H4V,SGN-B7H4V,SGNB 7H4V
Synonyms	VTCN1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	N/A
Host Species	Homo sapiens
IgG type	IgG1
Reactivity	Human
Target	B7-H4
Uniprot ID	Q7Z7D3
Description	PE-conjugated Anti-B7-H4(SGN-B7H4V biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

