

PRODUCT INFORMATION

Clone ID	DM181
Target	BAFF
Synonyms	TNFSF13B;BAFF;BLYS;CD257;DTL;TALL1;THANK;TNFSF20;ZTNF4;TALL-1
Host Species	Rabbit
Description	PE-conjugated Anti-BAFF antibody(DM181); Rabbit mAb
Delivery	3-4 weeks
Uniprot ID	Q9Y275
IgG type	Rabbit IgG
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	The protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This cytokine is a ligand for receptors TNFRSF13B:TACI; TNFRSF17:BCMA; and TNFRSF13C:BAFFR. This cytokine is expressed in B cell lineage cells; and acts as a potent B cell activator. It has been also shown to play an important role in the proliferation and differentiation of B cells. Alternatively spliced transcript variants encoding distinct isoforms have been identified.
Usage	Research use only
Conjugate	PE-conjugated

