

PRODUCT INFORMATION

Common Name	BIVV-020;BIVV020;BIVV020;SAR-445088;SAR445088;SAR445088;TNT-020;TNT020;TNT 020
Synonyms	EDSPD2
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	C15
Uniprot ID	P09871
Description	PE-conjugated Anti-C15(riliprubart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

