

PRODUCT INFORMATION

Common Name	IPH-5401, NN-8210, NNC-0215-0384
Synonyms	C5a-R;C5aR;CD88
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG1
Reactivity	Human
Target	C5AR1
Uniprot ID	P21730
Description	PE-conjugated Anti-C5AR1(avdoralimab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

