

PRODUCT INFORMATION

Common Name	LFG-316;LFG316;LFG 316;NOV-4;NOV 4;NOV4
Synonyms	C5D;C5a;C5b;ECLZB;CPAMD4
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	C5
Uniprot ID	P01031
Description	PE-conjugated Anti-C5(tesidolumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

