

PRODUCT INFORMATION

Common Name	KIND-016;KIND016;KIND 016;KIND-16;KIND16;KIND 16
Synonyms	Interleukin-31;IL-31
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb (Canine)
IgG type	IgG2
Reactivity	Canine
Target	Canine IL31
Uniprot ID	C7G0W1
Description	PE-conjugated Anti-Canine IL31(tirnovetmab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

