

PRODUCT INFORMATION

Common Name	KIND-016;KIND016;KIND 016;KIND-16;KIND16;KIND 16
Conjugate	PE-conjugated
Synonyms	Interleukin-31;IL-31
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	N/A
IgG type	IgG2
Reactivity	Canine
Target	Canine IL31
Uniprot ID	C7G0W1
Description	PE-conjugated Anti-Canine IL31(tirnovetmab biosimilar) mAb
Delivery	Under Development
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

