

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | 4F12 |
| Synonyms | PD1, PDCD1, CD279, SLEB2, hPD-1, hPD-l, hSLE1 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | N/A |
| Host Species | Chimeric Rabbit |
| IgG type | IgG |
| Reactivity | Canine |
| Target | PD-1 |
| Uniprot ID | NP_001301026.1 |
| Description | PE-conjugated Anti-Canine PD-1(INTERVET 4F12) mAb |
| Delivery | 3-4 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

