

PRODUCT INFORMATION

Common Name	IMGN-632;IMGN632;IMGN 632
Synonyms	IL3R;CD123;IL3RX;IL3RY;IL3RAY;hIL-3Ra;IL-3R-alpha
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD123
Uniprot ID	P26951
Description	PE-conjugated Anti-CD123(pivekimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

