

PRODUCT INFORMATION

Common Name	IMT009;IMT-009;IMT 009
Synonyms	NKR;CD161;CLEC5B;NKR-P1;NKR-P1A;hNKR-P1A
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD161
Uniprot ID	Q12918
Description	PE-conjugated Anti-CD161(rezorstobart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

