

PRODUCT INFORMATION

Common Name	ICP-B04;ICP B04;MOR 208;MOR208;MOR-208;XENP-5574;XENP5574;XENP 5574;XmAb-5574;XmAb5574;XmAb 5574
Synonyms	CD19;B4;CVID3;MGC12802
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD19
Uniprot ID	P15391
Description	PE-conjugated Anti-CD19(tafasitamab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

