

PRODUCT INFORMATION

Common Name	LY3454738;LY-3454738;LY 3454738
Synonyms	OX2R;MOX2R;CD200R;HCRTR2
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	CD200R1
Uniprot ID	Q8TD46
Description	PE-conjugated Anti-CD200R1(ucenprubart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

