

PRODUCT INFORMATION

Common Name	SCT-400;SCT400;SCT 400
Synonyms	B1;S7;Bp35;CD20;FMC7;CVID5;LEU-16;MS4A1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD20
Uniprot ID	P11836
Description	PE-conjugated Anti-CD20(ripertamab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

