

PRODUCT INFORMATION

Clone ID	DMC280
Target	CD21
Synonyms	C3DR; CD21; CR; CVID7; SLEB9
Host Species	Rabbit
Description	PE-conjugated Anti-CD21 antibody(DMC280); IgG1 Chimeric mAb
Delivery	3-4 weeks
Uniprot ID	P20023
IgG type	Rabbit/Human Fc chimeric IgG1
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	This gene encodes a membrane protein; which functions as a receptor for Epstein-Barr virus (EBV) binding on B and T lymphocytes. Genetic variations in this gene are associated with susceptibility to systemic lupus erythematosus type 9 (SLEB9). Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Usage	Research use only
Conjugate	PE-conjugated

