

PRODUCT INFORMATION

Common Name	SGN-35, cAC10-Val-Cit-MMAE
Synonyms	TNFRSF8;CD30;DIS166E;Ki-1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric
IgG type	IgG1
Reactivity	Human
Target	CD30
Uniprot ID	P28908
Description	PE-conjugated Anti-CD30 (brentuximab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

